



Discover the
INTIMITY MODE
of Doctor Tecar Therapy

Versatility, precision and safety.



AN EFFECTIVE SOLUTION FOR PELVIC PAIN RELIEF

INDICATIONS

- Pelvic perineal pain
- Lacerations
- Edemas, haematomas
- Adhesions
- Vaginismus
- Vulvodynia
- Urinary incontinence
- Trigger points
- Dyspareunias
- Vaginal dryness

*Doctor Tecar Therapy from Mectronic, thanks to its **frequency swing** and **polymodal emission**, is particularly effective in postnatal recovery and the **treatment of various conditions** such as vulvodynia, dyspareunia, incontinence, prolapse, vaginismus and pelvic floor dyssynergia.*

High-frequency energy stimulates cellular activity, **increasing metabolic activity and accelerating tissue regeneration**, improving vascular and lymphatic system involvement.



Childbirth and **menopause** can cause **many pathologies** that compromise women's quality of life.

Disorders such as **urinary incontinence, vaginal dryness, dyspareunia** and many others can be **improved by treatment** with *Doctor Tecar Therapy*.

INNOVATIVE TREATMENTS THANKS TO THE INTIMITY MODULE

APPLICATION BENEFITS OF THE INTIMITY MODULE

- Reduces **pain** and **urinary incontinence**.
- Improves **vaginal dryness** and **libido**.
- Tones the **perineum** and **prevents prolapse**.
- Promotes **post-partum healing**.
- Improves the **quality of sexual life**.

Interchangeable and **sterilisable heads** provide comfort and safety for the patient.



Equipped with a **sterilisable applicator**, it is effective in the management of **injuries** and promotes faster recovery.



Safe and effective, it revolutionises the treatment of pelvic floor disorders, improving quality of life.

DISCOVER THE MECTRONIC DOCTOR TECAR LINE

DOCTOR TECAR PLUS

7 Frequencies

300kHz | 448kHz | 500kHz |
750kHz | 1MHz | 12.MHz | 1.5MHz

Biotransfer system
with Biotransfer race

AVX automatic
sequential mode

220W
of effective power

3 emission modes
Continuous, Adjustable Pulsed and
Dimmer

Swing mode



DOCTOR TECAR SMART

2 Frequencies

448kHz | 750kHz

Innovative
Biotransfer System

AVX automatic
sequential mode

70 W
of effective power

2 emission modes
Continuous, Adjustable Pulsed

