



BIOSONYC

Biosonyc is a professional device for ultrasound therapy, characterized by the simplicity of use and the double emission frequency (1/3 MHz). Both probe in endowment and device can issue the suitable frequencies for all types of superficial and deep applications (both therapeutic and aesthetic).

Biosonyc contains pre-set programmes with the parameters for the most common applications and works by power mains; besides it is endowed with a comfortable small bag for the transport and the home therapy. Thanks to the Modul Line system it is possible to combine Biosonyc with up to 2 other devices among Magnetotherapy, Electrotherapy and Ionophoresis devices.



Programmes / Applications

Bone Atrophy	Tendonitis	Local Scleroderma
Bone Callus	Cicatrisation	Paravertebral Scleroderma
Osteitis and Periostitis	Oedematous Localised Cellulite	Warts
Delayed Ossification	Flabby Localised Cellulite	Bruises
Arthritis / Arthrosis	Massive Localised Cellulite	Tenosynovitis
Joint Stiffness	Oedematous Widespread Cellulite	Epicondylitis
Sprains	Flabby Widespread Cellulite	Periarthritis
Lumbago	Massive Widespread Cellulite	Bursitis
Myalgia	Localised Lymphatic Drainage	
Spasms	Widespread Lymphatic Drainage	
Muscular Tears	Foruncle	
Contractures	Granuloma	

Technical Specifications

Power Supply	AC 230V 50/60Hz
No. Output Channels	1
No. of Programmes	40 pre-set and 5 free
Frequency (Hz)	1/3 MHz
Emission	Continuous, modulated (1/2)
Power	0.1 - 3 W/cm ²
Dimensions and weight	22*17*17cm, 1,2 kg only device, 2 kg with accessories

Standard Kit

Bag
 1/3 MHz ultrasound probe Ø 45 mm
 Gel for ultrasound therapy 250g
 Feeding cable
 User's manual with photos of the applications