



BIOPHYSIO

Professional electrostimulator for physiotherapy; it issues all kind of electrotherapy waves for analgesic and anti-inflammatory applications, post-injury recovery of healthy and denervated muscles, urogynecological therapy. The software includes 50 pre-set programmes; it works by power means. Thanks to the Modul Line system, it is possible to combine Biophysio with up to 2 other devices for Magnetotherapy and Ultrasound therapy.



Programmes / Applications

Cervical Pain	Trophic Effect for Injuries with Oedemas
Periarthritis	Muscular Warm-up
Lumbago	Recovery from Fatigue
Epicondylitis	Recovery from Muscular Fatigue
Carpal Tunnel Syndrome	Post-injury Functional Rehabilitation
Knee Pain	Toning
Strain	Analgesic Urogynecological Applications
Tendonitis	Hyperhidrosis Applications
Acute Pains	Vasodilation and Hyperemising Effect
Chronic Pains	Stimulation of Denervated Muscles
Anti-inflammatory and Decontracting Effect	

Technical Specifications

Constant Current Generator	yes, 1	Standard Kit
Power Supply	power mains AC 220V 50/60Hz	Bag
Wave Forms	Tens, diadynamic, 2-pole interferential, rectangular biphasic, Kotz, ionophoresis, Faradic, rectangular monophasic, triangular, exponential, uro	2 electrode sponges 60x80mm
No. Output Channels	2 independent	2 electrode sponges 80x120mm
No. of Programmes	135 pre-set, 30 free	2 elastic bands
Intensity (mA)	0-120	4 self-adhesive plug electrodes 50x50
Frequency (Hz)	1-160	Feeding cable
Impulse Duration	50 µs - 1 sec	User's manual with photos of the applications
Dimensions and weight	22*17*17cm, 1,2 kg only device, 2 kg with accessories	