

FARMATEK L-1[®]

Farmatek L-1 is the new device for **CPM (Continuous Passive Motion)** for **rehabilitation of ankle and knee joint**.

Speed traction, angle, extension and time can be set according to patient's needs.

Endowed with a **remote controller**, it allows to safely stop and restart the therapy and to reverse traction direction.

Memory function: saves each selected parameter for the following session.

Large, easy-to-read-**LCD screen**.

Simple, user-friendly **control panel**.



APPLICATIONS

- ✓ Artropaties
- ✓ Artroscopies
- ✓ Pseudoarthrosis
- ✓ Endoprotheses
- ✓ Reconstruction ligaments
- ✓ Prophylaxis for thrombosis
- ✓ Fractures
- ✓ Arthritis
- ✓ Periarthritis

BENEFITS

- ✓ It accelerates the process of recovery
- ✓ It improves the mobility of the joints
- ✓ It reduces the hospital stay
- ✓ Good results of rehabilitation in shorter time and in a more comfortable way
- ✓ No side effects

TECHNICAL FEATURES

Moving distance of platform for knee: 0-410 mm

Moving distance of platform for ankle: 0-170 mm

Adjusting range of supporting frame: 300-620 (mm)

Extension range for knee: -5 - 120°

Extension range for ankle: -30°-0-45°

Platform's traction: 180-840 (N)

Supplier: power supply

50Hz / 220V / 80VA

Weight: 21 KG.

Dimensions: 820x280x300 mm

Treatment time: 0-250 (min.)

Moving speed platform: 0-10.5 (mm/s)